

DASH 6™

Improving patient safety in neuraxial connections.

ISO 80369-6 NRFit™ neuraxial connectors

NRFit™ & DASH 6™

- ISO 80369 establishes unique designs for different circuit applications.
- Each application is identified by its ISO 80369 suffix number.
- Connectors compliant with ISO 80369-6 are specific to neuraxial devices.
- The term “NRFit™” refers to neuraxial devices that are ISO 80369-6 compliant.
- Major manufacturers are using the colour yellow to identify NRFit™ devices.
- NRFit™ devices from GBUK Healthcare are exclusively branded “DASH 6™”.

 **GBUK
Healthcare**
Part of the GBUK Group 

Surety™

Surety™ neuraxial devices had previously been introduced to provide an alternative to the use of universal luer-type systems for neuraxial applications.

Neuraxial

NRFit™

DASH 6™ NRFit™ - Product Codes

Lock Syringes

Direct Code	NHSSC Code	Product Description	Unit of Sale
15011010	FWC475	1ml Lock Syringe	100pcs./box
15011030	FWC479	3ml Lock Syringe	100pcs./box
15011050	FWC483	5ml Lock Syringe	100pcs./box
15011100	FWC485	10ml Lock Syringe	100pcs./box
15011200	FWC487	20ml Lock Syringe	50pcs./box
15011600	FWC488	60ml Lock Syringe	25pcs./box



Accessories

Direct Code	NHSSC Code	Product Description	Unit of Sale
15241001	FWC491	3-Way Stopcock	50pcs./box
15111001	FWF149	Needle Hub Cap	100pcs./box
15041001	FWC489	Syringe Cap	100pcs./box
15171001	FWF150	Syringe to Syringe Adapter	100pcs./box
15231022	FWC490	0.22µ Disc Filter	100pcs./box
15151001	FTR2379	Blood Patch Kit	5pcs./box



ISO 80369-6:2018

is the international standard requirements for small-bore syringes intended to be used for neuraxial applications.

Neuraxial Applications

involve the use of medication to administer medication for wound infiltration anaesthesia, regional anaesthesia procedures or remove cerebro-spinal fluid for diagnostic purposes.

NRFit™

is the dedicated connector for regional block devices and syringes with luer connectors.

Dash 6™

is the exclusive brand name for the connectors from GBUK H.



Slip Syringes

Direct Code	NHSSC Code	Product Description	Unit of Sale
15021010	FWC477	1ml Slip Syringe	100pcs./box
15021030	FWC478	3ml Slip Syringe	100pcs./box
15021050	FWC481	5ml Slip Syringe	100pcs./box
15021100	FWC484	10ml Slip Syringe	100pcs./box
15021200	FWC486	20ml Slip Syringe	50pcs./box
15021101	Direct Only	10ml Slip LOR Device	50pcs./box



10ml Slip LOR Device



016

ard that specifies
ore connectors
connections

ications

al devices intended
s to neuraxial sites,
hesis delivery, and other
cedures or to monitor
al fluid for therapeutic

tor for neuraxial and
nd is not compatible

me for the NRFit™
Healthcare.

Direct Code	NHSSC Code	Product Description	Unit of Sale
15311022	FWC493	0.22µ Epidural Flat Filter	100pcs./box



15281100	FWF151	5µ, 100mm Drawing Up Filter Straw	100pcs./box
----------	--------	-----------------------------------	-------------



15491050	FTR2344	18g x 50mm, Blunt Drawing Up Needle without Filter	100pcs./box
----------	---------	--	-------------

15291050	FTR2343	18g x 50mm, Blunt Drawing Up Needle with 5µ Filter	100pcs./box
----------	---------	--	-------------




15461010	FWL106	10cm Extension Set with Clamp	25pcs./box
----------	--------	-------------------------------	------------



15461125	FWL107	125cm Extension Set with Clamp	25pcs./box
----------	--------	--------------------------------	------------

15161001	FTR2380	Ommaya Reservoir Kit	5pcs./box
----------	---------	----------------------	-----------



Surety™ to DASH 6™ Conversion Chart **NRFit™**

Surety™ Code	Product Description	DASH 6™ NRFit™ Code
05011010	1ml Lock Syringe	15011010
05011030	3ml Lock Syringe	15011030
05011050	5ml Lock Syringe	15011050
05011100	10ml Lock Syringe	15011100
05011200	20ml Lock Syringe	15011200
05011600	60ml Lock Syringe	15011600
05021010	1ml Slip Syringe	15021010
05021030	3ml Slip Syringe	15021030
05021050	5ml Slip Syringe	15021050
05021100	10ml Slip Syringe	15021100
05021200	20ml Slip Syringe	15021200
05021101	LOR Device 10ml Slip	15021101
05041001	Syringe Cap	15041001
05111001	Needle Hub Cap	15111001
05241001	3-Way Stopcock	15241001
05281100	100mm, 5µ Drawing Up Filter Straw	15281100
05291050	18g x 50mm, Blunt Drawing Up Needle with 5µ Filter	15291050
N/A	18g x 50mm, Blunt Drawing Up Needle without Filter	15491050
05311022	0.22µ Epidural Flat Filter	15311022
05231022	0.22µ Bacterial Disc Filter	15231022
05171001	Syringe to Syringe Adapter	15171001
05461010	10cm, Extension Set with Clamp	15461010
05461125	125cm, Extension Set with Clamp	15461125
05151000	Blood Patch Kit	15151001
05161001	Ommaya Reservoir Kit	15161001



GBUK Group award-winning specialist companies are expert in the provision of enteral feeding, critical care and patient moving and handling products and technologies. GBUK Group companies are suppliers to every UK NHS Trust and to a rapidly expanding Export customer base.

gbukgroup.com



gbukgroup.com/nrfit



GBUK Healthcare
Woodland House
Blackwood Hall Business Park
North Duffield, Selby
North Yorkshire YO8 5DD

T +44 (0) 1757 288 587
E info@gbukhealthcare.com